



Patient name _____ Date _____

Device OrthoApnea NOA OrthoApnea Classic

Patient measurements

Example

Maximum retrusion -4 mm

Maximum protrusion 6 mm

Advancement range |-4|+6|= 10 mm

Starting point (SP) 60% de 10 = 6 -> - 4 + 6 = +2 mm

Dental midline deviation in occlusion
 To the right To the left / None

Dental midline deviation in protrusion
 To the right To the left / None

OrthoApnea NOA device configuration settings Optional

Vertical dimension Minimum distance between posterior cusps: 3mm. As per registration Minimum possible

Advancement sequence type Measurements indicated with respect to SP. They can be indicated in "mm" or in "%". Standard -1 SP +1 +2 Customized SP Additional splints

Device opening limit Ranging from 0 - 12mm. By default: 7mm.


Laterality Ranging from 0 to 5mm. By default: 4mm.

Buccal bands Upper Lower

Settings Anterior opening Hooks

Design of the splints Scallop Mixed scallop

Observations



Mark those teeth where you wish to alleviate retention